

NEW



OXY'PHARM®

NOCOSPRAY®

new generation

- One device to treat volumes from 10 to 1 000 m³
- Setting by touchpad
- Strengthened design



www.oxypharm.net



TECHNICAL SPECS

Cycle management	Touch pad NEW
Mist particle size	Around 5 µm
Decontamination ratio	From 1 ml/m ³
Room size/volume treatment	10 m ³ to 1000 m ³
Dimensions (mm)	342 H x 494 W x 300 D NEW
Net weight	6 kg
Power supply	220-240 V / 50 Hz - 1000 W 230-240 V / 60 Hz - 1000 W 115 V / 50Hz-60Hz - 1000 W

Use biocides safely. Always read the label and product information before use.



WEIGHS ONLY 6KG



Setting by touch pad

NEW

treats volumes from 10 to 1000 m³ with a single device

NEW

To be used with OXY'PHARM products : Nocolyse and Oxyxy range.



ECONOMICAL, FAST ACTION

Automated system where 1ml of product is sufficient for treating 1m³.



LIGHT WEIGHT AND COMPACT

The Nocospray is device for processing rooms with volumes of up to 1000 m³ and weighs only 6Kg.



FOR DISINFECTION: FULL ACTIVITY SPECTRUM

Bactericidal, virucidal, fungicidal and sporicidal, as per the current standards.



NORMS

Validated according to several norms, including NF T 72281 standard (2014)



SAFE

Non-allergenic products without risk to the environment or to humans.



FOR INSECTICIDE TREATMENT : BEST QUALITY/PRICE RATIO

Proven and approved effectiveness against bed bugs with costs among the lowest on the market.



OXY'PHARM®

917 rue Marcel Paul • 94500 Champigny-sur-Marne • France
T : +33 1 48 82 58 29 • F : +33 1 48 82 46 13 • commercial@oxypharm.net

www.oxypharm.net

